

Lastogrip®



Spec.-No.: **D 9.4013**
 Department: **CMO-DGM**
 Date: **2015-12-08**

1. General Product Description

Lastogrip® is an elasticated tubular bandage, exerts uniform, medium-strength pressure, quick and easy to apply without additional fixation, fit perfectly without slipping, may be removed as often as required and re-applied, kind to the skin due to a high content of cotton fibres; may be washed and sterilized.

Lastogrip® carries the CE mark according to EU directive 93/42/EEC for medical devices.

The product is classified as a class I medical device unsterile.

A conformity assessment has been performed for Lastogrip® and it has been shown to be in compliance with all applicable requirements of the above mentioned directive.

The safe use of Lastogrip® therefore, is ensured if the product is used in line with the intended purpose (see 2. Application/Indication).

2. Application / Indication

- Stronger compressing fixation bandages, e.g. in Thoracic and Abdominal surgery
- For post Mamma-OP bandages
- For wide abrasions to avoid oedema at the wound border respectively to relieve the wound
- As supporting and relieving bandage for damaged of the attitude and movement apparatus
- As supporting thermal bandage for arthrosis, and so on

3. Presentations

Lastogrip® natural white, rolls length approx. 10 m ± 0,5 m

Size	Width (cm) (unexpanded) approx.	Indication
Size A	ca. 4,20	for children's hand and arms
Size B	ca. 6,25	for small hands, arms and feet
Size C	ca. 6,75	for medium-sized hands, arms, feet and lower legs
Size D	ca. 7,50	for large hands, arms, feet and lower legs
Size E	ca. 8,75	for legs, large arms and small thighs
Size F	ca. 10,00	for legs and thighs
Size G	ca. 12,00	for large thighs
Size J	ca. 17,50	for small trunk
Size K	ca. 21,50	for small trunk
Size L	ca. 32,50	for large trunk
Size M	ca. 37,50	for very large trunk

Lastogrip®



Spec.-No.: **D 9.4013**
 Department: **CMO-DGM**
 Date: **2015-12-08**

Lastogrip® C skin-colored, rolls length approx. 10m ± 0,5 m

Size	Width (cm) (un-expanded) approx.	Indication
Size B	ca. 6,25	for small hands, arms and feet
Size C	ca. 6,75	for medium-sized hands, arms, feet and lower legs
Size D	ca. 7,50	for large hands, arms, feet and lower legs
Size E	ca. 8,75	for legs, large arms and small thighs
Size F	ca. 10,00	for legs and thighs
Size G	ca. 12,00	for large thighs

4. Product Characteristics

Material composition:

Cotton: approx. 82 %
 Elastodiene: approx. 10 % (contains latex)
 Polyamide: approx. 8 %

Right-right binding on round knitting machine.
 Elastic element at every 4th row.

5. Product requirements

Elongation (cross direction): > 250 %
 Sterilizability: Steam according to DIN EN ISO 17665-1 with validated process

6. Labelling

No. with 9-digit code:

Lot. - No.	5	999	01	00	1
Explanation	[5 = 2015]	[999 - 000]	Internal Code		
Definition	Year	Serial production order			

The expiry date is printed near the symbol ⌘.

Expiry date ⌘ Year Month
 Example 2020 12
 Shelf life 5 years

7. Packaging

Lastogrip® or Lastogrip® C rolls are individually packed into dispensers. The dispensers are packed into transport cartons, sealed with adhesive tape and stacked onto Euro pallets.

Lastogrip®



Spec.-No.: **D 9.4013**


Department: **CMO-DGM**

Date: 2015-12-08

Date: 12th December 2015

PAUL HARTMANN AG
Product Development CMO-DGM

i.A.

A handwritten signature in blue ink, appearing to read "C. Jäger".

Cathrin Jäger