		<h1>Technical Data Sheet</h1>
--	--	-------------------------------

<b>Foliodrape® Protect Plus, Hip Set II reinforced (Art. No. 938 772/1)</b>		Spec.-No.:	D 6.1001-01
	Department:	CMO-DOE	
	Date:	2015-11-02	

## 1. General Product Description

*Foliodrape® Protect Plus Hip Set II reinforced* for hip surgeries like femoral nailing, hip endoprosthesis, etc.

The *Foliodrape® Protect Plus Hip Set II reinforced* carries the CE-mark according to EU directive 93/42/EWG for medical devices and is classified as class I sterile medical device.

A conformity assessment has been performed for *Foliodrape® Protect Plus* and it has been shown to be in compliance with all applicable requirements of the above mentioned directive.

The safe use of *Foliodrape® Protect Plus Hip Set II reinforced*, therefore, is ensured if the product is used in line with the intended purpose.


## 2. Application / Indication

*Foliodrape® Protect Plus Hip Set II reinforced* is a disposable product for sterile draping of the patient and the theatre equipment.

## 3. Presentations

<b>Art. No. 938 772/1</b>	<b><i>Foliodrape® Protect Plus Hip Set II reinforced</i></b>
	sterile, 3 Sets per dispenser, 1 dispenser per transport carton
	Dimensions of dispenser: 545 x 374 x 234 mm (L x W x H)
	Dimensions of transport carton: 590 x 390 x 250 mm (L x W x H)

## 4. Product Characteristics

<p><b>Content:</b></p> 	1 reinforced table cover 140 x 190 cm 1 mayo stand cover (telescopic folded) 80 x 145 cm 1 self-adhesive drape 170 x 300 cm (self-adhesive 100 cm) 1 self-adhesive, reinforced slit sheet drape 280 x 245 cm (slit 10 x 95 cm, reinforced zone 150 x 160 cm) 1 approach drape 140 x 190 cm 1 drape 75 x 90 cm 1 stockinet 35 x 120 cm 4 cellulose towels 30 x 33 cm 2 adhesive tapes 10 x 50 cm <u>underlined</u> is the side with the adhesive zone on
<p><b>Sterilization:</b> <b>Self-adhesives:</b></p>	With ethylene oxide gas according to EN 556-1 and EN 11135-1 Skin-friendly synthetic rubber, latex free and colophony free

<b>Foliodrape® Protect Plus, Hip Set II reinforced (Art. No. 938 772/1)</b>		Spec.-No.:	D 6.1001-01
		Department:	CMO-DOE
		Date:	2015-11-02

5. Product Requirements

*Foliodrape® Protect Plus* corresponds to all requirements of European standard for surgical drapes EN 13795, for the highest performance level “High performance, critical area”.

For Details see Technical Data sheet of the base material (D 6.1005a, D 6.1005b), mayo stand cover (D 6.1005c), back table cover (D 6.1005d) and stockinet (D 6.5417a).

6. Labelling

Lot-No. with 9-digit code:

e.g.: 5023XXXXX

LOT	5	023	XXXXX
	year of production	serial production order	for internal use

Expiry date:

e.g. ⌚ 2020 – 01  
year month

Shelf life: 5 years

The product will be sterile during this period as long as the packaging is not damaged and the storage conditions are kept.

7. Packaging

*Foliodrape® Protect Plus Hip Set II reinforced* is sealed into sterile packaging (peeling pouch) according to EN ISO 11607 and EN 868. The peel pouch is packed into folding box and transport carton, sealed with adhesive tapes and stored on euro-pallet.

Date: 2015-11-02

PAUL HARTMANN AG  
Development OTM & Exam. Gloves (CMO-DOE)

i. A.

  
Dr. Evelyn Hartmann