

## MoliCare® Premium MEN PAD 3 Drops/green



Spec.-No.: D 9.0411

Department: COH-IDT

Date: 2018-04-20

### 1. General Product Description

Incontinence pad specially designed for men; with Dry Plus Distribution Layer and superabsorbent core; neutralises odour; complete nonwoven covering; active backsheet with textile comfort; novel allround Protect Cuff System; provides optimum security by highest absorption and storing capacity and guarantees protection against re-wetting; anatomic V – shape provides a comfortable fit.

Classification acc. EEC-Directive for non sterile medical devices class I

### 2. Application / Indication

MoliCare MEN Pads 3 Drops are absorbent liners for incontinence care for adults. Except MoliCare MEN PAD 2 Drops which is for male adults. The products are suitable for usage in case of bladder weakness with several severities.

### 3. Presentations

MoliCare® Premium MEN PAD 3 Drops

### 4. Product Characteristics

<b>Material:</b>	<b>MoliCare® Premium MEN PAD 3 Drops</b>
- liquid protection:	Nonwoven-film-laminate, air active
- absorbent core:	100 % kraft pulp, bleached elementary chlorine free, kraft pulp with superabsorbent polymer
- distributor edition:	Nonwoven, white
- adhesive:	Hotmelt adhesive
- nonwoven:	Polypropylene nonwoven, hydrophilic, white
- cuffs:	Polypropylene nonwoven hydrophobic, white
- elastics:	Elastics for standing gathers
- backing:	Silicon paper


## MoliCare® Premium MEN PAD 3 Drops/green



Spec.-No.: D 9.0411  
Department: COH-IDT  
Date: 2018-04-20

### 5. Product Drawing

MoliCare Premium MEN PAD 3 Drops



Dimensions [mm]	MoliCare Premium MEN PAD 3 Drops
total length:	$260 \pm 5\%$
total width:	$180 \pm 5\%$
Total length absorbent core:	$215 \pm 5\%$
total width absorbent core abdominal	$140 \pm 5\%$

## MoliCare® Premium MEN PAD 3 Drops/green



Spec.-No.: D 9.0411  
Department: COH-IDT  
Date: 2018-04-20

### 6. Product Requirements

Requirements:		MoliCare Premium MEN PAD 3 Drops
Total weight*:	[g]	15 ± 10%

\* ± 10% Tolerance corresponds with ± 3s

### 7. Labeling

Lot-No. with 6-digit code:

e.g.:

XXX	XX	X
day	year	for internal purposes only
1 - 365	20 <u>12</u>	1

### 8. Packaging

MoliCare® Premium MEN PAD 3 Drops are single wrapped and packed into polyethylene-bag; into transport carton acc. DIN, packed onto europallet.

April 20<sup>th</sup> 2018



- Product Development COH-IDT -